

Tips and Syringes



Black Micro FX
Black Mini
Black Mini Brush
Black Micro
Blue Micro
Blue Mini Dento-Infusor
Capillary
Chromaclone PVS Mixing
Intraoral Impression
Deliver-Eze Brush
Endo-Eze
Endo-Eze Irrigator
FX Flex
Inspiral Brush
Metal Dento-Infusor
Micro Capillary
Micro 20ga

Micro 20ga FX
NanoTip
NaviTip
NaviTip FX
SurgiTip
Mixing
Intraoral
White Mac
White Mini Laser
Skini Syringes
1.2, 3, and 5ml Syringes
Ultradent Syringe Covers
Luer Lock Caps
Syringe Organizers
Ultradent Micro Applicators
TriAway Adapter



Don't settle for wannabes.



Often imitated but never duplicated -

Ultradent's direct delivery tips streamline your procedures by eliminating extra steps and reducing waste. Backed by three decades of quality and reliability, Ultradent's tips offer a range of choices and features.

Check out our tips with **LOK-TITE®**

Luer Lock tips with Lok-Tite feature double threads that lock the tip into place for increased security, and wings for easy attachment and removal.



22ga

Black Micro® FX® Tips

- Accommodates various viscosities
- Flocked tip fans out to spread materials in a thin, uniform layer

Optimal with: PrepQuick and PrimaDry.

LOK-TITE®

#1357 - 100pk

#1434 - 500pk



Black Mini® Tips

- Dispenses large volumes
- Opaque plastic preserves the flow of light-cured materials

Optimal with: Ultra-Blend plus, Ultradent LC Block-Out, PermaFlo, PermaFlo DC, TriAway Adapter, UltraTemp, Opalescence Boost, OpalDam, and Opalescence Endo.

LOK-TITE®

#514 - 100pk

#1433 - 500pk



Black Mini® Brush Tips

- Provides precise, controlled delivery of aqueous materials
- Tight, adjustable brush fibers minimize bubbles
- Unique to Ultradent

Optimal with: Consepsis, Diamond Polish Mint, Peak SE, PermaFlo DC Primers, PermaQuick Primers, PermaSeal, Seek/Sable Seek, Ultradent Silane, and Universal Dentin Sealant.

LOK-TITE®

#1169 - 100pk

#193 - 20pk w/tipsocs

#1432 - 500pk

Length of fibers can be customized on brush tips by grasping bristles and moving them in or out of tip. Shortening bristles gives a stiffening effect; lengthening bristles gives a softening effect and also allows the fibers to fan out and cover larger areas.



Black Micro[®] Tips
Blue Micro[®] Tips

- Provides pinpoint precision
- Narrow cannulae accurately deliver materials

BLACK MICRO optimal with: Ultra-Blend plus, EtchArrest, OpalDam, UltraSeal XT plus, and Ultracare.

BLUE MICRO optimal with: Ultra-Etch.

#1085 - Black 100pk
#127 - Blue 100pk

#1435 - Black 500pk
#1436 - Blue 500pk



Blue Mini[®]
Dento-Infusor[®] Tips

- Offers the same tissue management benefits as the Metal-Dento Infusor
- Allows controlled flow of drop-sized quantities

Optimal with: Astringedent, Astringedent X, and ViscoStat Regular, Wintermint, and Clear.

LOK-TITE[®]

#1086 - 100pk

#1440 - 500pk w/tipsocs



Capillary Tips

- Navigates canals and substantially minimizes use of paper points
- Narrow, flexible taper accesses curved canals

Optimal with: Luer Vacuum Adapter for moisture removal from Endodontic canals.

LOK-TITE[®]

#341 - 0.014" 20pk
#186 - 0.019" 20pk

#3099 - 0.014" 50pk
#1425 - 0.019" 50pk



Chromaclone[®] PVS
Mixing Tips

- Enables direct delivery of impression materials
- Auto-mixing, disposable, and color-coded

Optimal with: Chromaclone PVS.

#2904 - Green 50pk
#2902 - Yellow 50pk

#2903 - Pink 50pk



IntraOral
Impression Tips

- Allows precise placement
- Attaches to PVS mixing tip

Optimal with: Chromaclone PVS.

#2906 - 50pk

→ Tips and Syringes



Deliver-Eze® Brush Tips

- Allows precise placement
- Only bendable metal cannulae brush tip that fits into post preparations as small as 0.6mm in diameter

Optimal with: PermaFlo DC Primers A and B.

#5125 - 50pk

#5126 - 100pk



Endo-Eze® Tips

- Works ideally for a variety of Endodontic procedures
- Comes with flexible, strong cannulae in a variety of sizes
- Bends easily

Optimal for: luting materials, and air/water delivery.

Optimal with: TriAway Adapter, PermaFlo DC, and other Ultradent Syringes.

20pk

100pk

Gauge - Diameter

#345

#1428

18ga - 0.049";1.25mm

#346

#1429

19ga - 0.042";1.06mm

#347

#1430

20ga - 0.035";0.90mm


#348

#1431

22ga - 0.028";0.70mm

#350

30ga - 0.010";0.25mm




Endo-Eze® Irrigator Tips

- Provides ideal reach without risk of expressing chemicals past the apex
- Comes with a flexible, blunt cannula with a unique, anti-obturating end
- 1.0" non-sterile and 1.25" sterile

Optimal with: ChlorCid, ChlorCid V, and Ultradent 5ml syringe.

#207 - 20pk

#202 - Sterile 20pk



FX® Flex Tips

- Provides visible, controlled delivery
- Flexible neck and brush tip spread materials evenly

Optimal with: Flor-Opal Varnish White.

LOK-TITE®

#6025 - 25pk

#6026 - 75pk



Inspiral® Brush Tips

- Delivers viscous or filled materials smoothly via an internal helical channel and ridge
- Tight, adjustable brush fibers minimize bubbles

Optimal with: Composite Wetting Resin, HiSite, Peak LC, PermaQuick Resin, PQ Amalgam, PQ1, Ultradent Porcelain Etch, Ultra-Etch and UltraSeal XT plus.

LOK-TITE®

#123 - 100pk

#1033-I - 500pk



19ga

Metal Dento-Infusor® Tips

- Places hemostatic agents precisely and effectively removes superficial coagulum
- Blunt, bent cannula with padded brush enables gentle pressure on the sulcus
- Ultradent's first tip, the "MDI" remains paramount for successful tissue management

Optimal with: Astringedent, Astringedent X, and ViscoStat Regular, Wintermint, and Clear.

LOK-TITE®

#2558 - 20pk

#2559 - 100pk

#2560 - 500pk



5mm

10mm

Micro Capillary™ Tips

- Bright colors are easily identified against soft tissues
- The world's smallest molded tip (inside diameter: 0.008")

Optimal with: Ultracare, perio materials, endodontics, and Luer Vacuum Adapter.

LOK-TITE®

#1120 - 5mm 20pk

#1802 - 5mm 100pk

#1121 - 10mm 20pk



Micro 20 Gauge Tips

- Large gauge cannula enables consistent flow
- Standard flowable composite delivery tip

Optimal with: PermaFlo, PermaFlo Purple, and OpalDam.

LOK-TITE®

#1252 - 100pk

#1437 - 500pk



Micro 20 Gauge FX® Tips

- Standard flowable delivery tip
- Flocked fibers fan out to evenly spread material over the entire surface

Optimal with: Opalescence Boost.

LOK-TITE®

#1405 - 100pk

#1406 - 500pk

→ Tips and Syringes



NanoTip® Tips

- Ideal for micro dentistry
- Small, strong, and rigid 30 gauge cannula

Optimal with: UltraSeal XT Plus, Ultra-Etch, PermaFlo, etc., microscopic restorative dentistry, and delivery of Vit-I-escence Colors for esthetic characterizations.

#1240 - 20pk

#1426 - 100pk



NaviTip® Tips

- Provides controlled delivery to the apex
- Flexible, stainless steel cannula easily navigates curved canals
- Available in 17, 21, 25, and 27mm lengths and 29 or 30 gauge cannulae

Optimal with: EndoREZ, Ultradent EDTA 18% Solution, ChlorCid, ChlorCid V, Consepsis, Consepsis V, Ultradent Citric Acid 20% Solution, File-Eze, and UltraCal XS.


	29ga-20pk	30ga-20pk	29ga-50pk	30ga-50pk
17mm, White	#5112	#1249	#1378	#1421
21mm, Yellow	#5113	#1349	#1374	#1422
25mm, Blue	#5114	#1250	#1376	#1423
27mm, Green	#5115	#1354	#1377	#1424
Assorted	#5116	#1351	#1379	#3319

NaviTip® 31ga with Double Sideport Irrigator Tip

- Double sideports deliver irrigants safely, minimizing the possibility of chemicals being expressed past the apex
- The world's smallest cannula navigates the most intricate canal spaces

Optimal with: ChlorCid, EDTA Solution, Citric Acid 20% Solution, and Consepsis.

Note: 31ga should not be used with viscous solutions.



31ga-20pk
#5121 - 21mm
#5123 - 27mm

31ga-50pk
#5122 - 21mm
#5124 - 27mm



NaviTip® FX® Tips

- One-of-a-kind brush cleans, scrubs, and irrigates simultaneously
- Available in 17 or 25mm lengths, with a rigid 30 gauge cannula

Optimal with: Ultradent Citric Acid 20% Solution.

LOK-TITE®

#1452 - 17mm 20pk
#1454 - 25mm 20pk

Listed as an "EXCELLENT" product by a prominent independent research institute.¹

¹ "Can't Live Without" Clinical Research Associates Newsletter, Page 4, January, 2005



SurgiTip™ Tips

- Ideal for delicate surgeries
- Large diameter tip opening

Optimal with: Luer Vacuum Adapter for small periodontic or endodontic surgeries and controlled suction of Opalescence Boost.

LOK-TITE®

#1248 - 20pk



Mixing Tips

- Mixes and delivers in one action

Optimal with: Ultradent Dual-Barrel syringes - UltraTemp, EndoREZ, and PermaFlo DC.

#5920 - Mixing 20pk

#5922 - Intraoral 20pk



Intraoral Tips

- Allows precise placement
- Attaches to Dual-Barrel Mixing Tip

Optimal with: PermaFlo DC.



White Mac™ Tips and White Mini™ Laser Tips

- Dispenses large volumes
- All-plastic delivery tips

Optimal with: OraSeal Caulking, OraSeal Putty, Consepsis Scrub, and Opalustre.

#1361 - 100pk
White Mac

#1480 - 500pk
White Mac

#1247 - 20pk
White Laser



Skini Syringes

- Generates pressure in the syringe with low force to the plunger, resulting in more precise apical delivery
- Low waste
- Content: 0.5ml

Optimal with: EndoREZ and PermaFlo DC.

#1680 - 20pk

#1681 - 50pk



1.2ml Syringes

- Snug; Luer Lock threads prevent tips from popping off
- Designed for intraoral use (shapes and smoothes with rounded edges and contours)

#124 - 20pk

#157 - 100pk



3ml Impression Syringes

- Syringe barrel flanges positioned for optimum control/leverage

Optimal with: Any Ultradent tips and elastomeric impression materials.

#600 - 10pk



5ml Syringes

- Syringe barrel flanges positioned for optimum control/leverage

Optimal with: Irrigants for in-office or dentist-supervised procedures, as well as Capillary tip, ChlorCid, Citric Acid EDTA, and Consepsis.

#201 - 10pk



Ultradent® Syringe Covers

- Provides an easy, reliable barrier
- Ensures asepsis of syringe during clean-up

Optimal with: All 1.2ml syringes

#249 - 250pk

Tips and Syringes



Luer Lock Caps

- Winged, polypropylene, plastic Luer Lock caps
- Use to seal syringes loaded in the office

Optimal with: All Ultradent plastic syringes.

#205 - 20pk



Full Size Syringe Organizer

- Holds 30 syringes
- Sliding, 12-compartment tray
- Made of injection-molded polycarbonate

#379 - Full Size Organizer



Half Size Syringe Organizer

- Holds 14 syringes
- Made of clear acrylic

#382 - Half Size Organizer

Ultradent® Micro Applicators

- Come with standard size tips
- Bend easily
- Apply small amounts accurately

The micro applicator is comprised of small non-linting, non-absorbent fibers that suspend solutions, creating a non-dripping sphere.



#1106 - Micro Applicators 400pk - 200 each color



TriAway™ Adapter

Use for controlled administration of water and/or air to depths of tiny cavity preparations, such as minimally invasive operative-type preparations or endodontic preparations. Do NOT use in open canals.



1. Wash preparations.



2. Rinse and dry etchants.



3. Thin and air dry bonding resins.

#1375 - TriAway Disposable Irrigators 20pk

